

STARKEY PRODUCT LINE



STARKEY PRODUCT LINE

If you are having technical difficulties, stay logged in and call Audiology Online at 1-800-753-2160.





STARKEY PRODUCT LINE

- This session is worth 1 hour / .1 CEU credit
- Must stay logged in for the entire session
- Must successfully complete a short quiz









NOISE MANAGEMENT





ACOUSTIC SCENE ANALYZER SOLVING THE PROBLEM OF NOISE BY PRIORITIZING SPEECH INVISION DIRECTIONALITY Improves SNR Reduces listening effort Fast-acting system Provides listener comfort Reduces listening effort Reduces listening effort

VOICE IQ2

- Significantly improves speech clarity even in the noisiest environments
- Reduces listening effort and cognitive fatigue
- Instantly applies variable noise adaptation in all channels between each pause in speech
- Provides up to 20dB of noise reduction
- Available in all devices in the 110, 90 and 70 technology levels in:





INVISION DIRECTIONALITY

- · Industry leading directivity
- · Varying low frequency compensation
- Dynamic switching between Omni and Directional
- Prioritizes speech understanding when patients need it most



Available in all product families and technology levels



AUDIOSCAPES





Available in all product families and technology levels



BINAURAL SPATIAL MAPPING

- Boosts the performance of our Acoustic Scene Analyzer: provides best SNR possible!
- Delivers maximum listening comfort and provides a more natural, 3D-like sound quality
- Minimizes patient's perception of disparity and switching







IQ BOOST

- Available with SurfLink Mobile OR the Advanced SurfLink Remote
- Allows patients to quickly enable aggressive noise control performance
- Sound quality and speech is preserved
- Available in all wireless technology levels in:









DIRECTIONALITY PLUS

- Provides maximum adaptation for Voice iQ² when directionality is engaged (up to 20dB)
- Adjusts the input level required for maximum activation from 85 dB SPL down to 80 dB SPL
- Available in all wireless technology levels in these families:

 3series XIDO








FEEDBACK CANCELLATION



PUREWAVE FEEDBACK ELIMINATOR 🔥



- Proven up to 28 dB of added stable gain in an open fit
- Proven stability in every hearing aid fit
- Available in all devices in the 110, 90, 70 and 30 technology levels in:








Starkey. Hearing Is Our Concern

FREQUENCY LOWERING



SPECTRAL IQ

- · Intelligently identifies high-frequency speech cues, then replicates them in lower frequencies
- Features the broadest active bandwidth of any lowering technology
- · Causes less distortion than other approaches due to preservation of spectral cues and harmonic relationships







SPECTRAL IQ

WWW.SPECTRALIQ.COM



Available in all products and technology 3 series Xino Lijseries X series tour ignite





WIRELESS SOLUTIONS





SURFLINK ACCESSORIES 100% wireless Proprietary 900 mHz IRIS wireless platform Compatible with ALL Starkey wireless devices

SURFLINK MOBILE



Wireless. Seamless. Mobile. One Device does it all.

- Cell phone streaming with JUST TALK!
- · Media signal streaming
- Remote microphones (ALD function)
- · Hearing aid remote control





INVISIBLE SOLUTIONS



INVISIBLE OPTIONS





STARKEY PEDIATRICS





STARKEY PEDIATRICS

- Full product line
- Bight colors
- Dynamic FM
- Pediatric Auto Path
- 1 year of free earmold remakes
- Exclusive partnership with Nickelodeon
- Starkey Pediatrics Mobile App!






PROFESSIONAL & PATIENT TOOLS







HEARING LOSS SIMULATOR



Inspire Software



LIFESTYLE SOLUTIONS







Inspire SoftwareApple App



LIFESTYLE SOLUTIONS



Starkey.

NEW! EXPERIENCE MANAGER





Available in all Inspire-based technology families



SWAP FIT

- · Transfer patient settings from previous product to a new one
- · Helps provide familiar starting point



Available in all Inspire-based technology families



SPEECH MAPPING

2



- Speech Mapping screen Access from any of these programming
- - **Quick Fit**
- Fine Tuning
- Compression
- Spectral iQ Memories
- Expert Assistant

Available in all Inspire-based technology families





- · Patient partners in the adjustment process
- · Research shows that outcomes are improved when patients are personally involved in the decision making process
- · Inspire or iPad applications
- · Easy to use and understand
- Available in all devices in the 110, 90 and 70 technology levels in:





INDEPENDENT SPEECH OPTIMIZATION

- Allows for independent control of soft, moderate and loud level speech
- Adjust for specific patient reports and let the compression take care of itself!
- Offers a more personalized fitting approach to optimize the listening experience
- Available in all technology levels in:







SELF LEARNING



Gradually and automatically learns patient volume control preferences in each active memory.

Normal memory: 12 hour learning window Other memories: 2 hour learning window Available in all devices in:







SYNCHRONIZED USER ADJUSTMENTS



NON-ACOUSTIC TELEPHONE & STREAMING MEMORIES

 Makes audibility easier when using telephones, telecoils, or telecoil loops and offers higher fidelity sound when enjoying streamed audio. · Streaming memory is available in all wireless products in: 3series Xino wisenes |

tour Ignite

· Telephone memories enhanced in all families



TELECOIL LOOP MEMORY

- · Adds a 12 dB boost to the telecoil setting
- · Can be implemented binaurally
- · Appropriate for patients who use loop systems in public places



Available in all product families and technology levels supporting a telecoil



WIRELESS DEMOS

Features

- Flagship technology level
- Omnidirectional Microphone
- Wireless
- 115/50 matrix
- Automatic Telephone Response (ATR)
- 312 battery











THERE'S AN APP FOR THAT!























XPERIA











HEAR COACH

- Listening games that challenge both cognitive and
- Track progress over time and unlock levels
- · Lower return for credit
- · Improve hearing instrument benefit
- · Increase patient delight
- · Lower perceived Hearing Handicap

FREE to download!



STARKEY ON DEMAND

T² Remote T² On Demand Audiology On Demand











LIFESTYLES AND LEVELS

Lifestyle	Technology Level	Channels/Bands
VIBRANT	i110/110	16
ACT I V E	i90/90	12
SOCIAL	i70/70	8
CALM	i30/30	6
BASIC	i20/20	4

Wireless Technology: 'i' precedes the technology level














Highlights

- Rocker Switch for volume and memory
- Advanced case design to improve directivity
- · Snap Fit receivers
- Independent Speech Optimization
- Cell phone and media connectivity (wireless)
- Improved non-acoustic input Streaming (wireless) and Telephone (wireless + wired) memories
- Directionality Plus (wireless)
- · Automatic DAI functionality








RIC DESIGN & DURABILITY

- Snap Fit receiver Secure six pin connector Oval shape more resistant to twisting
- Wrap-around mic and switch cover
- Improved battery door design
 Case design promotes airflow over the external surface while creating a natural handling area
- Improved RIC cable cosmetics







3 series

NEW BTE 13!

- · A small and stylish wireless BTE with the longevity of a 13 battery
- Features 3 Series' patient-driven design enhancements
- Wireless and wired devices in three technology levels and price points: 110, 90 & 70.









Standard Features

Voice iQ² InVision Directionality AudioScapes Spectral iQ Pure Wave Feedback Eliminator

Telecoil / Autocoil

ISO

SoundPoint Speech Mapping Swap Fit Self Learning T2 on Demand Audiology on Demand HydraShield²



Wireless	Wired
RIC 312*	RIC 312
BTE 13	BTE 312
BTE 312	



3 series



Wireless Advantages

Directionalilty Plus Binaural Spatial Mapping iQ Boost Synchronized User Adjustments Streaming Memory

Cable and Dongle Free Programming Media & Cell phone Connectivity Intelligent Media Mobility





3 series











Economy Wireless and Wired

Tour i30/30: 6 channels/bands Tour i20/20: 4 channels/bands







RIC 312* RIC 312* BTE 312



Standard Features

Select Noise Adaptation InVision Directionality AudioScapes Spectral iQ Pure Wave Feedback Eliminator

Telecoil / Autocoil Speech Mapping Swap Fit T2 on Demand Audiology on Demand



BTE 312 BTE 13





Wireless Advantages

- Synchronized User Adjustments
- Cable and Dongle Free Programming
- Media & Cell phone connectivity • Intelligent Media Mobility
- · Streaming Memory























Highlights

- Incredibly small, wireless micro RIC 312
- Fully configurable, ergonomic switch
- Easy volume and memory adjustments
- Cell phone and media connectivity Same great feature set at 3 Series Wireless!









Highlights

- Micro RIC 10
- 5 Technology Levels: 16, 12, 8, 6 and 4 channels
- Programmable Capacitive Switch VC or memory
- Induction Coil
- HydraShield²
- Standard and Bright Colors
- International Markets













Highlights

- · Micro RIC 10
- Tinnitus Therapy + Amplification in one device!
- · Configurable Memories
 - Tinnitus Therapy only
 - Tinnitus Therapy + Amplification
 - Amplification only





Highlights

- Active Lifestyle * 90 Technology Level
- Best fit protocol
- Allows for flexible adjustments in 16 bands
- · Amplification independent of tinnitus stimulus
- · Adjustable tinnitus stimulus modulation rate
- · Flexible device-level control of volume and memory
- Can adhere to any fitting philosophy, including TRT and masking protocols













SoundPoint Tinnitus

- · Allows patients to identify the most preferable stimulus
- · Ensures satisfaction
- · Inspire and iPad Applications
- Easy to use!








































Compatible with all SurfLink Accessories









EDUCATION AND TRAINING



EDUCATIONAL OPPORTUNITIES

StarkeyPro.com StarkeyEvidence.com Hearing-Aid.com QuickTIPS
Technology Summits
Regional Experience Summits

SpectraliQ.com Live Online Courses
Inspire Online Recorded Online Courses
You Tube Podcasts Personal Training



Call your Starkey Sales Representative for more information about any of these opportunities!



